

## Realize Language Starter's Guide

### Creating a customized target vocabulary

Being able to use a variety of parts of speech is a basic skill in sentence building. Realize Language provides simple tools to help you support your student to learn this skill. In Analysis you have seen the default vocabulary list in Word Groups. This is a list created by the well-known Speech Pathologist Gail van Tatenhove. The list provides is a fairly extensive goal for student to work towards. You can use it to generate your own goals for specific parts of speech. In this example we will use adjectives, but a goal can be created for any part of speech.

- Step 1: Select adjectives, verbs and prepositions. You will be able to observe your student's performance in each group with this default list.
- Step 2: Note any words that you want to track for your client
- Step 3. Go to manage goals and create a new list
- Step 4: Enter your customized list and select store.
- Step 5: Return to Word groups.

In the coming weeks and months, you will be able to follow your students increasing use of these words. Realize Language also allows you to share student's data with other account holders. In this way, you will be able to support working with your student.

