

Realize Language Starter's Guide

Activity Four: Comparing word lists

Select the **Words** tab and then the **List** section at the bottom of the page.

The Word list section allows you to examine the words used by your student and compare words that were correctly spelled or attempted and those that are found in the extensive Realize Language dictionary (e.g. *cat*) with those that are not (e.g. *Antonio*).

Explore

- Select the four buttons KNOWN, UNKNOWN, PRESTORED and SPELLED and observe how the list changes.
- Check the number next to the word. This tells you how frequently the word was used.
- Set the Lists to UNKNOWN and PRESTORED. Now you are looking at any pre-stored words used that are unknown to the database, for instance people's names.
- Set the List to SPELLED and KNOWN and then to UNKNOWN. You will now see any words that were spelled letter by letter.



Questions

Should the names of any friends or family members be stored into the device?

Can you spot words your student was trying to spell?

Try This

- Press the export button to create two text lists. One of all the pre-stored words that have been used and a second of all the known spelled words. Print the lists and mark any spelled words which are actually stored in the device.
- Begin a list of words that you think you might store into the device.
- Begin a list of words that persons modeling on the device should be able to easily find.



Known	Unknown	Pre-stored	Spelled	Export
154 i	6 two	3 asked	1 next	1 maybe
61 to	6 eat	3 listening	1 from	1 ahead
38 drink	6 away	3 orange	1 paper	1 me
37 want	6 ask	3 boy	1 goldsmith	1 god
37 yogurt	6 bear	3 multiplicity	1 i'm	1 spaghetti
37 no	6 office	3 am	1 same	1 world
35 at	6 monday	3 food	1 number	1 church
32 yellow	6 out	2 where	1 children	1 original
28 star	5 thursday	2 everyone	1 bad	1 over
26 yes	5 kiwi	2 around	1 gold	1 activate
26 after	5 sad	2 them	1 rain	1 first
26 the	5 don't	2 mixed	1 bed	1 reason
24 music	5 know	2 cinnamon	1 colors	1 sunday
22 on	5 goldfish	2 cherry	1 clown	1 very
21 card	5 procedures	2 blast	1 toy	1 so
21 and	5 mountain	2 ice	1 what	1 answers
21 right	5 favorite	2 cream	1 visits	1 doctor
19 home	5 class	2 story	1 bring	1 today
19 now	5 boat	2 because	1 xerox	1 snake
18 my	5 fall	2 up	1 warm	1 yell
18 strawberry	5 thank	2 his	1 look	1 other
18 play	5 was	2 different	1 ready	1 chickens
17 toolbox	5 are	2 green	1 october	1 afraid
17 go	5 saw	2 see	1 points	1 top
16 listen	5 she	2 car	1 dream	1 lunch
15 minute	5 watch	2 best	1 dreaming	1 boost
15 monsters	5 goodbye	2 subtract	1 something	1 bonus
14 is	4 have	2 water	1 has	1 material

Known	Unknown	Pre-stored	Spelled	Export
11 antonio	1 pgg	1 cac	1 tokyo	1 hidhid
11 inc	1 za	1 pointv	1 toon	1 mnet
10 vivaldi	1 scatch	1 pushingis	1 ity	1 movfishing
5 boyer	1 scrach	1 20seteps	1 vvivaldi	1 bszip
3 mrsrms	1 goki	1 er	1 hfklixbear	1 movwe
3 tex	1 3zac's	1 crbab	1 cbacdb	1 dtffinss
2 mrs	1 cacbcc	1 cocah	1 owl	1 peedy
2 york	1 cba-	1 yyork	1 oeowl	1 ionation
2 aa	1 xpop	1 ne	1 apeothe	1 nerig
2 scarch	1 il	1 lefaf	1 toldchic	1 di
2 li	1 foagain	1 ppelab	1 mismiss	1 beivh
2 paul	1 herar	1 hhiiii	1 smiss	1 monser
2 mike	1 antonnio	1 aantonio	1 trl	1 tawfanmat
2 disney	1 vi	1 acaa	1 homrework	1 kstream
2 ms	1 antantonio	1 sno	1 wot	1 jane
2 nigel	1 lis	1 snones	1 whot	1 utnto